

PON2 Rabbit mAb


ABZ2104

Product Overview


Name	PON2 Rabbit mAb
Accession(Primary)	Q15165
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 39 kDa; Observed MW: 40-42 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; IP-1:50-1:100

Target information (Q15165)

Gene Name	PON2
About Gene Symbol	PON2 (Paraoxonase 2) is a protein phosphatase involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of PON2 in K562, HeLa lysates using PON2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using PON2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.