


PKC mu/PKD Rabbit mAb

ABZ2106

Product Overview

Name	PKC mu/PKD Rabbit mAb
Accession(Primary)	Q15139
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 102 kDa; Observed MW: 102 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q15139)

Gene Name PRKD1**About Gene Symbol**

Western blot analysis of PKC mu/PKD in C6, HeLa lysates using PKD1 antibody.

