


## PDK2 Rabbit mAb

ABZ2107

### Product Overview

|                             |                                              |
|-----------------------------|----------------------------------------------|
| <b>Name</b>                 | PDK2 Rabbit mAb                              |
| <b>Accession(Primary)</b>   | Q15119                                       |
| <b>Host Species</b>         | Rabbit                                       |
| <b>Clonality</b>            | Monoclonal                                   |
| <b>Isotype</b>              | IgG                                          |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                            |
| <b>Calculated MW</b>        | Calculated MW: 46 kDa; Observed MW: 46 kDa   |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IHC-P-1:200-1:1000; IP-1:20 |

### Target information (Q15119)

**Gene Name** PDK2**About Gene Symbol**

Western blot analysis of PDK2 in rat Brain lysates using PDK2 antibody.

