

PDK1 Rabbit mAb


ABZ2108

Product Overview

Name	PDK1 Rabbit mAb
Accession(Primary)	Q15118
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 49 kDa; Observed MW: 49 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:10-1:100


Target information (Q15118)

Gene Name	PDK1
About Gene Symbol	PDK1 (Pyruvate Dehydrogenase Kinase 1) is a metabolic dehydrogenase involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright© 2024 Zen-bioscience Limited.

Western blot analysis of PDK1 in LNCaP lysates
using PDK1 antibody.



Western blot analysis of PDK1 in LNCaP lysates
using PDK1 antibody.