


## RAGE Rabbit mAb

ABZ2109

### Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | RAGE Rabbit mAb                            |
| <b>Accession(Primary)</b>   | Q15109                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Mouse, Rat                                 |
| <b>Calculated MW</b>        | Calculated MW: 43 kDa; Observed MW: 55 kDa |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000                           |

### Target information (Q15109)

**Gene Name** AGER**About Gene Symbol**

Western blot analysis of RAGE in mouse lung  
lysates using RAGE antibody.

