

p40-phox Rabbit mAb


ABZ2110

Product Overview


Name	p40-phox Rabbit mAb
Accession(Primary)	Q15080
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 39 kDa; Observed MW: 39 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200; FC-1:20-1:100; IP-1:20-1:50

Target information (Q15080)


Gene Name	NCF4
About Gene Symbol	



Western blot analysis of p40phox in HL-60, U251 lysates using NCF4 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using p40phox antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of p4phox (green)
in hela using p4phox antibody, and DAPI(blue)