

## P2Y6 Rabbit mAb


ABZ2111

### Product Overview


|                             |                                                                          |
|-----------------------------|--------------------------------------------------------------------------|
| <b>Name</b>                 | P2Y6 Rabbit mAb                                                          |
| <b>Accession(Primary)</b>   | Q15077                                                                   |
| <b>Host Species</b>         | Rabbit                                                                   |
| <b>Clonality</b>            | Monoclonal                                                               |
| <b>Isotype</b>              | IgG                                                                      |
| <b>Species Reactivity</b>   | Human                                                                    |
| <b>Calculated MW</b>        | Calculated MW: 36 kDa; Observed MW: 36 kDa                               |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:100-1:200;<br>ICC/IF-1:100-1:200; IP-1:20-1:50 |

### Target information (Q15077)


|                          |       |
|--------------------------|-------|
| <b>Gene Name</b>         | P2RY6 |
| <b>About Gene Symbol</b> |       |



Immunocytochemistry analysis of P2Y6 (green) in A549 using P2Y6 antibody, and DAPI (blue).



Western blot analysis of P2Y6 in K562, HeLa lysates using P2Y6 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer tissue using P2Y6 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.