

Periostin Rabbit mAb


ABZ2113

Product Overview


Name	Periostin Rabbit mAb
Accession(Primary)	Q15063
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 93 kDa; Observed MW: 93 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:2000; IP-1:20-1:50

Target information (Q15063)

Gene Name	POSTN
About Gene Symbol	



Western blot analysis of Periostin in C6 lysates using Periostin antibody.



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using Periostin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.