

ACAP2 Rabbit mAb


ABZ2114

Product Overview

Name	ACAP2 Rabbit mAb
Accession(Primary)	Q15057
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 88 kDa; Observed MW: 88 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:100-1:200

Target information (Q15057)

Gene Name	ACAP2
About Gene Symbol	ACAP2 (ArfGAP With Coiled-Coil, Ankyrin Repeat And PH Domains 2) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of ACAP2 in HeLa, CHO-K1, C6, Jurkat lysates using ACAP2 antibody.