

NOLC1 Rabbit mAb


ABZ2117

Product Overview


Name	NOLC1 Rabbit mAb
Accession(Primary)	Q14978
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 74 kDa; Observed MW: 110 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:200-1:500; FC-1:200-1:500

Target information (Q14978)


Gene Name	NOLC1
About Gene Symbol	



Western blot analysis of NOLC1 in Jurkat lysates using NOLC1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using NOLC1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of NOLC1 (green)
in HeLa using NOLC1 antibody, and DAPI (blue)