

STAT4 Rabbit mAb


ABZ2118

Product Overview


Name	STAT4 Rabbit mAb
Accession(Primary)	Q14765
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 86 kDa; Observed MW: 86 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:500-1:2000; ICC/IF-1:20-1:50; FC-1:50-1:100; IP-1:20-1:50

Target information (Q14765)

Gene Name	STAT4
About Gene Symbol	



Western blot analysis of STAT4 in HeLa, A549 lysates using STAT4 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using STAT4 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.