

MVP Rabbit mAb


ABZ2119

Product Overview

Name	MVP Rabbit mAb
Accession(Primary)	Q14764
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 99 kDa; Observed MW: 110 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:50-1:200; ICC/IF-1:100-1:500; FC-1:100-1:500

Target information (Q14764)

Gene Name	MVP
About Gene Symbol	MVP (Major Vault Protein) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of MVP in U87-MG lysates using MVP antibody.