

Lamin B Receptor Rabbit mAb


ABZ2120

Product Overview


Name	Lamin B Receptor Rabbit mAb
Accession(Primary)	Q14739
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Rat
Calculated MW	Calculated MW: 71 kDa; Observed MW: 71 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:50-1:100; ICC/IF-1:200-1:500; FC-1:200-1:1000

Target information (Q14739)

Gene Name	LBR
About Gene Symbol	



Western blot analysis of Lamin B Receptor in K562, HeLa lysates using Lamin B Receptor antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Lamin B Receptor antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.