

USP10 Rabbit mAb


ABZ2122

Product Overview


Name	USP10 Rabbit mAb
Accession(Primary)	Q14694
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 87 kDa; Observed MW: 110 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:200-1:500; IP-1:10-1:100

Target information (Q14694)


Gene Name	USP10
About Gene Symbol	



Western blot analysis of USP10 in HeLa lysates using USP10 antibody.



Immunocytochemistry analysis of USP10(green) in HeLa using USP10 antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using USP10 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.