

SMC1A Rabbit mAb

ABZ2124


Product Overview

Name	SMC1A Rabbit mAb
Accession(Primary)	Q14683
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 143 kDa; Observed MW: 143 kDa
Recommended Dilution	WB-1:2000-1:5000; IHC-P-1:200-1:500

Target information (Q14683)

Gene Name SMC1A

About Gene Symbol



Western blot analysis of SMC1A in K562, rat Brain, C6, 3T3, HeLa lysates using SMC1A antibody.

