

IRF3 Rabbit mAb


ABZ2126

Product Overview


Name	IRF3 Rabbit mAb
Accession(Primary)	Q14653
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 47 kDa; Observed MW: 47 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:100-1:200; FC-1:20; IP-1:20

Target information (Q14653)

Gene Name	IRF3
About Gene Symbol	IRF3 (Interferon Regulatory Factor 3) is an interferon regulatory transcription factor involved in JAK-STAT/cytokine signaling and innate immune signaling. It is an important research target in immunology, inflammation, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of IRF3 in HeLa, A549, HL-60 lysates using IRF3 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using IRF3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.