

IL-11RA Rabbit mAb


ABZ2127

Product Overview


Name	IL-11RA Rabbit mAb
Accession(Primary)	Q14626
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 45 kDa; Observed MW: 45 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-1:100-1:200; IF-1:50-1:100; FC-1:50-1:100

Target information (Q14626)

Gene Name	IL11RA
About Gene Symbol	



Western blot analysis of IL-11RA in K562 lysates using IL-11 Receptor alpha antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using IL-11RA antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.