

Desmocollin 3 Rabbit mAb


ABZ2129

Product Overview


Name	Desmocollin 3 Rabbit mAb
Accession(Primary)	Q14574
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 100 kDa; Observed MW: 100 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200; ICC/IF-1:50-1:100

Target information (Q14574)

Gene Name	DSC3
About Gene Symbol	



Immunocytochemistry analysis of Desmocollin 3 (green) in HeLa using Desmocollin 3 antibody, and DAPI (blue).



Western blot analysis of Desmocollin 3 in K562, HeLa lysates using Desmocollin 3 antibody.