

MCM6 Rabbit mAb


ABZ2130

Product Overview


Name	MCM6 Rabbit mAb
Accession(Primary)	Q14566
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 93 kDa; Observed MW: 93 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:500-1:2000; ICC/IF-1:200-1:500; FC-1:50-1:100; IP-1:20-1:50

Target information (Q14566)


Gene Name	MCM6
About Gene Symbol	MCM6 (Minichromosome Maintenance Complex Component 6) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechani



Western blot analysis of MCM6 in Jurkat lysates using MCM6 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using MCM6 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of MCM6 (green)
in HeLa using MCM6 antibody, and DAPI (blue)