

ENT2 Rabbit mAb


ABZ2131

Product Overview


Name	ENT2 Rabbit mAb
Accession(Primary)	Q14542
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 50 kDa; Observed MW: 60 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200; FC-1:20

Target information (Q14542)


Gene Name	SLC29A2
About Gene Symbol	



Western blot analysis of ENT2 in K562 lysates using ENT2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using ENT2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of ENT2 (green) in
Hela using ENT2 antibody, and DAPI(blue)