

Metabotropic Glutamate Receptor 2 Rabbit mAb


ABZ2133

Product Overview


Name	Metabotropic Glutamate Receptor 2 Rabbit mAb
Accession(Primary)	Q14416
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 96 kDa; Observed MW: 99 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:2000; ICC/IF-1:200-1:500; FC-1:50-1:100

Target information (Q14416)

Gene Name	GRM2
About Gene Symbol	



Immunohistochemistry analysis of paraffin-embedded mouse testis using Metabotropic Glutamate Receptor 2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of Metabotropic Glutamate Receptor 2 in K562, HeLa lysates using Metabotropic Glutamate Receptor 2 antibody.