

GNA13 Rabbit mAb

ABZ2134


Product Overview


Name	GNA13 Rabbit mAb
Accession(Primary)	Q14344
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 44 kDa; Observed MW: 44 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000

Target information (Q14344)

Gene Name GNA13

About Gene Symbol





Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using G protein alpha 13 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.