


PTK2B Rabbit mAb

ABZ2136


Product Overview

Name	PTK2B Rabbit mAb
Accession(Primary)	Q14289
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 116 kDa; Observed MW: 116 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:100-1:200

Target information (Q14289)

Gene Name PTK2B**About Gene Symbol**

Western blot analysis of PYK2 in A549, HL-60, U2OS lysates using PYK2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using PYK2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.