

## TRIM25 Rabbit mAb


ABZ2137

### Product Overview


|                             |                                                                                           |
|-----------------------------|-------------------------------------------------------------------------------------------|
| <b>Name</b>                 | TRIM25 Rabbit mAb                                                                         |
| <b>Accession(Primary)</b>   | Q14258                                                                                    |
| <b>Host Species</b>         | Rabbit                                                                                    |
| <b>Clonality</b>            | Monoclonal                                                                                |
| <b>Isotype</b>              | IgG                                                                                       |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                                                         |
| <b>Calculated MW</b>        | Calculated MW: 71 kDa; Observed MW: 71 kDa                                                |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:100-1:200;<br>ICC/IF-1:200-1:1000; FC-1:50-1:100; IP-1:20-1:100 |

### Target information (Q14258)

|                          |                                                                                                                                                                                                                                                                 |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | TRIM25                                                                                                                                                                                                                                                          |
| <b>About Gene Symbol</b> | TRIM25 (Tripartite Motif Containing 25) is a ubiquitin E3 ligase involved in ubiquitin-proteasome pathway. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and i |



Western blot analysis of TRIM25/EFP in K562, rat Brain, C6, HeLa lysates using TRIM25 antibody.



Immunocytochemistry analysis of TRIM25/EFP (green) in MCF-7 using TRIM25/EFP antibody, and DAPI (blue)