

Cortactin Rabbit mAb


ABZ2138

Product Overview


Name	Cortactin Rabbit mAb
Accession(Primary)	Q14247
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 62 kDa; Observed MW: 80-85 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:500-1:1000; FC-1:20; IP-1:20-1:50

Target information (Q14247)


Gene Name	CTTN
About Gene Symbol	



Western blot analysis of Cortactin in K562, rat Brain, C6, 3T3, Hela lysates using Cortactin antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using Cortactin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of
Cortactin(green) in Hela using Cortactin
antibody,and DAPI(blue)