

Keap1 Rabbit mAb


ABZ2139

Product Overview


Name	Keap1 Rabbit mAb
Accession(Primary)	Q14145
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 70 kDa; Observed MW: 60-64 kDa
Recommended Dilution	WB-1:2000-1:10000; IHC-1:50-1:100

Target information (Q14145)


Gene Name	KEAP1
About Gene Symbol	KEAP1 (Kelch Like ECH Associated Protein 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Keap1 in Raw264.7 lysates using Keap1 antibody.



Western blot analysis of Keap1 in HeLa, A549, C2C12 lysates using Keap1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using Keap1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.