

IL-18 Rabbit mAb


ABZ2140

Product Overview


Name	IL-18 Rabbit mAb
Accession(Primary)	Q14116
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 22 kDa; Observed MW: 17 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:1000; FC-1:20-1:50

Target information (Q14116)

Gene Name	IL18
About Gene Symbol	



Western blot analysis of IL-18 in HeLa lysates using IL-18 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using IL-18 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.