

Protein CASP Rabbit mAb


ABZ2143

Product Overview


Name	Protein CASP Rabbit mAb
Accession(Primary)	Q13948
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 77 kDa; Observed MW: 77 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:100-1:200; FC-1:50-1:100; IP-1:20-1:50

Target information (Q13948)


Gene Name	CUX1
About Gene Symbol	



Immunocytochemistry analysis of Protein CASP (green) in A549 using Protein CASP antibody, and DAPI(blue).



Western blot analysis of Protein CASP in 3T3, Raw264.7 lysates using Protein CASP antibody.



Western blot analysis of Protein CASP in HeLa, A549, HL-60, U2OS, C6 lysates using Protein CASP antibody