

NEAS Rabbit mAb


ABZ2144

Product Overview


Name	NEAS Rabbit mAb
Accession(Primary)	Q13813
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 285 kDa; Observed MW: 285 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q13813)

Gene Name	SPTAN1
About Gene Symbol	SPTAN1 (Spectrin Alpha, Non-Erythrocytic 1) is a structural cytoskeletal protein involved in cytoskeletal organization. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of NEAS in HeLa, Hacat lysates using NEAS antibody.



Western blot analysis of NEAS in mouse brain, rat brain lysates using NEAS antibody.