

Integrin alpha 9 Rabbit mAb

ABZ2145


Product Overview

Name	Integrin alpha 9 Rabbit mAb
Accession(Primary)	Q13797
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 114 kDa; Observed MW: 140 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200

Target information (Q13797)


Gene Name ITGA9

About Gene Symbol



Western blot analysis of Integrin alpha 9 in A549 lysates using Integrin alpha 9 antibody.

Copyright© 2024 Zen-bioscience Limited.



Western blot analysis of Integrin alpha 9 in HepG2 lysates using Integrin alpha 9 antibody.