

# Activin Receptor Type IIB Rabbit mAb

ABZ2148


## Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | Activin Receptor Type IIB Rabbit mAb       |
| <b>Accession(Primary)</b>   | Q13705                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                          |
| <b>Calculated MW</b>        | Calculated MW: 58 kDa; Observed MW: 58 kDa |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IP-1:20-1:50             |

## Target information (Q13705)

**Gene Name** ACVR2B

**About Gene Symbol**



Western blot analysis of Activin Receptor Type IIB/ACVR2B in CHO-K1, mouse Brain lysates using Activin Receptor Type IIB antibody.

