

## FHL1 Rabbit mAb


ABZ2150

### Product Overview


|                             |                                                                       |
|-----------------------------|-----------------------------------------------------------------------|
| <b>Name</b>                 | FHL1 Rabbit mAb                                                       |
| <b>Accession(Primary)</b>   | Q13642                                                                |
| <b>Host Species</b>         | Rabbit                                                                |
| <b>Clonality</b>            | Monoclonal                                                            |
| <b>Isotype</b>              | IgG                                                                   |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                                     |
| <b>Calculated MW</b>        | Calculated MW: 36 kDa; Observed MW: 36 kDa                            |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IHC-P-1:200-1:500;<br>ICC/IF-1:20-1:50; IP-1:20-1:50 |

### Target information (Q13642)

|                          |      |
|--------------------------|------|
| <b>Gene Name</b>         | FHL1 |
| <b>About Gene Symbol</b> |      |



Western blot analysis of FHL1 in HeLa, U251 lysates using FHL1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using FHL1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.