

## SLC12A1/NKCC2 Rabbit mAb

ABZ2151

### Product Overview

|                             |                                              |
|-----------------------------|----------------------------------------------|
| <b>Name</b>                 | SLC12A1/NKCC2 Rabbit mAb                     |
| <b>Accession(Primary)</b>   | Q13621                                       |
| <b>Host Species</b>         | Rabbit                                       |
| <b>Clonality</b>            | Monoclonal                                   |
| <b>Isotype</b>              | IgG                                          |
| <b>Species Reactivity</b>   | Human                                        |
| <b>Calculated MW</b>        | Calculated MW: 121 kDa; Observed MW: 150 kDa |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; FC-1:100-1:500             |

### Target information (Q13621)

**Gene Name** SLC12A1**About Gene Symbol**

