

Cullin 4A/4B Rabbit mAb


ABZ2152

Product Overview


Name	Cullin 4A/4B Rabbit mAb
Accession(Primary)	Q13620
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat, Hamster
Calculated MW	Calculated MW: 104 kDa; Observed MW: 88,104 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200

Target information (Q13620)

Gene Name	CUL4B
About Gene Symbol	



Western blot analysis of Cullin 4A/CUL4A+Cullin 4B/CUL4B in Jurkat, C6, CHO-K1 lysates using Cullin 4A/4B antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Cullin 4A/B antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.