


Cullin 3 Rabbit mAb

ABZ2154


Product Overview

Name	Cullin 3 Rabbit mAb
Accession(Primary)	Q13618
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 89 kDa; Observed MW: 89 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:20-1:100

Target information (Q13618)

Gene Name CUL3**About Gene Symbol**

Western blot analysis of Cullin 3/CUL3 in PC-12, rat spleen lysates using Cullin 3 antibody.



Western blot analysis of Cullin 3/CUL3 in K562, 3T3, HeLa lysates using Cullin 3/CUL3 antibody.