


Cullin 2 Rabbit mAb

ABZ2155

Product Overview

Name	Cullin 2 Rabbit mAb
Accession(Primary)	Q13617
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 87 kDa; Observed MW: 87 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:50-1:100; IP-1:10-1:100

Target information (Q13617)

Gene Name CUL2**About Gene Symbol**

Western blot analysis of Cullin 2/CUL2 in K562, C6, 3T3, HeLa lysates using Cullin 2 antibody.

