

IQGAP2 Rabbit mAb


ABZ2160

Product Overview


Name	IQGAP2 Rabbit mAb
Accession(Primary)	Q13576
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 181 kDa; Observed MW: 181 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100

Target information (Q13576)

Gene Name	IQGAP2
About Gene Symbol	IQGAP2 (IQ Motif Containing GTPase Activating Protein 2) is a GTPase-activating protein involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of IQGAP2 in A549, HL-60 lysates using IQGAP2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using IQGAP2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.