

IRF5 Rabbit mAb


ABZ2161

Product Overview

Name	IRF5 Rabbit mAb
Accession(Primary)	Q13568
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 56 kDa; Observed MW: 56 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:500; ICC/IF-1:500-1:1000; FC-1:20-1:50; IP-1:100-1:200


Target information (Q13568)

Gene Name	IRF5
About Gene Symbol	



Copyright© 2024 Zen-bioscience Limited.

Western blot analysis of IRF5 in Raw264.7 lysates using IRF5 antibody.



Western blot analysis of IRF5 in THP-1 lysates using IRF5 antibody.