

CaMKII delta Rabbit mAb


ABZ2163

Product Overview

Name	CaMKII delta Rabbit mAb
Accession(Primary)	Q13557
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 56 kDa; Observed MW: 56 kDa
Recommended Dilution	WB-1:500-1:1000; HC-P:1:50-1:100; IHC-1:50-1:100

Target information (Q13557)

Gene Name	CAMK2D
About Gene Symbol	CAMK2D (Calcium/Calmodulin Dependent Protein Kinase II Delta) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is wi



Western blot analysis of CaMKII delta in rat Brain, C6, HeLa lysates using CaMKII delta antibody.