

CaMKII Rabbit mAb


ABZ2164

Product Overview

Name	CaMKII Rabbit mAb
Accession(Primary)	Q13554
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 73 kDa; Observed MW: 54 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:100-1:200

Target information (Q13554)

Gene Name	CAMK2B
About Gene Symbol	



Western blot analysis of CaMKII in mouse brain
lysates using CaMKII beta antibody.