


RIP Rabbit mAb

ABZ2166

Product Overview

Name	RIP Rabbit mAb
Accession(Primary)	Q13546
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 76 kDa; Observed MW: 76 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q13546)

Gene Name RIPK1**About Gene Symbol**

Western blot analysis of RIP in K562, HeLa lysates using RIP antibody.

