

Pin1 Rabbit mAb


ABZ2167

Product Overview


Name	Pin1 Rabbit mAb
Accession(Primary)	Q13526
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 18 kDa; Observed MW: 18 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:100; FC-1:200-1:500; IP-1:20-1:50; IHC-1:50-1:100

Target information (Q13526)


Gene Name	PIN1
About Gene Symbol	




Immunocytochemistry analysis of PIN1 (green) in 293 using PIN1 antibody, and DAPI (blue).



Western blot analysis of Pin1 in mouse spleen lysates using Pin1 antibody.



Western blot analysis of Pin1 in Jurkat, rat Brain, C6, CHO-K1, Hela lysates using Pin1 antibody



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using Pin1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.