

cIAP1 Rabbit mAb


ABZ2172

Product Overview


Name	cIAP1 Rabbit mAb
Accession(Primary)	Q13490
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 70 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:50-1:100

Target information (Q13490)

Gene Name	BIRC2
About Gene Symbol	



Western blot analysis of cIAP1 in K562 lysates using cIAP1 antibody.



Immunocytochemistry analysis of cIAP1(green) in Hela using cIAP1 antibody, and DAPI(blue)