

Syntrophin alpha 1 Rabbit mAb


ABZ2176

Product Overview


Name	Syntrophin alpha 1 Rabbit mAb
Accession(Primary)	Q13424
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 54 kDa; Observed MW: 54 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:100-1:200; FC-1:100-1:200

Target information (Q13424)


Gene Name	SNTA1
About Gene Symbol	



Immunocytochemistry analysis of Syntrophin alpha 1 (green) in U87-MG using Syntrophin alpha 1 antibody, and DAPI (blue).



Immunohistochemistry analysis of paraffin-embedded mouse colon using Syntrophin alpha 1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of Syntrophin alpha 1 in
K562 lysates using Syntrophin alpha 1 antibody