

Ikaros Rabbit mAb


ABZ2177

Product Overview


Name	Ikaros Rabbit mAb
Accession(Primary)	Q13422
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 58 kDa; Observed MW: 50-70 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; FC-1:50-1:100

Target information (Q13422)

Gene Name	IKZF1
About Gene Symbol	



Western blot analysis of Ikaros in CHO-K1, Jurkat lysates using Ikaros antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Ikaros antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.