

ILK Rabbit mAb


ABZ2178

Product Overview


Name	ILK Rabbit mAb
Accession(Primary)	Q13418
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 51 kDa; Observed MW: 51 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200; IP-1:20

Target information (Q13418)

Gene Name	ILK
About Gene Symbol	ILK (Integrin Linked Kinase) is a serine/threonine protein kinase involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Integrin linked ILK in mouse lung lysates using Integrin Linked ILK antibody.



Western blot analysis of Integrin linked ILK in C6, HeLa lysates using Integrin linked ILK antibody.