

CTBP1 Rabbit mAb


ABZ2179

Product Overview


Name	CTBP1 Rabbit mAb
Accession(Primary)	Q13363
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 48 kDa; Observed MW: 48 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-1:100-1:200; IF-1:100-1:200; FC-1:200-1:2000

Target information (Q13363)


Gene Name	CTBP1
About Gene Symbol	



Immunocytochemistry analysis of CTBP1 (green)
in HEPG2 using CTBP1 antibody, and DAPI (blue).



Western blot analysis of CtBP1 in A549, MCF-7,
U251 lysates using CtBP1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using CtBP1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.