

ATM Rabbit mAb


ABZ2180

Product Overview


Name	ATM Rabbit mAb
Accession(Primary)	Q13315
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 351 kDa; Observed MW: 351 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:200-1:500; FC-1:20

Target information (Q13315)

Gene Name	ATM
About Gene Symbol	ATM (ATM Serine/Threonine Kinase) is a serine/threonine protein kinase involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunocytochemistry analysis of ATM (green) in HT-1080 using ATM antibody, and DAPI (blue).



Western blot analysis of ATM in HeLa, 293 lysates using ATM antibody.