

KAP1 Rabbit mAb


ABZ2183

Product Overview


Name	KAP1 Rabbit mAb
Accession(Primary)	Q13263
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 89 kDa; Observed MW: 100 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:100-1:200

Target information (Q13263)


Gene Name	TRIM28
About Gene Symbol	



Western blot analysis of KAP1 in K562, C6, 3T3, HeLa lysates using KAP1 antibody.



Immunocytochemistry analysis of KAP1 (green) in HeLa using KAP1 antibody, and DAPI (blue)



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using KAP1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.