

NME3 Rabbit mAb


ABZ2184

Product Overview


Name	NME3 Rabbit mAb
Accession(Primary)	Q13232
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 19 kDa; Observed MW: 19 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; IP-1:50-1:100

Target information (Q13232)

Gene Name	NME3
About Gene Symbol	



Western blot analysis of NME3 in A549, HL-60 lysates using NME3 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using NME3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.